

Baňka Erlenmeyerova úzkohrdlá

Kónická baňka s orientační stupnicí

- Bílá orientační stupnice a popisovací plocha
- Borosilikátové sklo



BORO 3.3						PYREX®					
Objem (ml)	Průměr baňky (mm)	ID hrdla (mm)	Výška (mm)	Balení (ks)	Katalogové číslo	Objem (ml)	Průměr baňky (mm)	ID hrdla (mm)	Výška (mm)	Balení (ks)	Katalogové číslo
5	25	12	52	1		5	25	12	52	1	S113001
10	30	12	60	1		10	30	12	60	1	S113004
25	42	22	75	1	B103006	25	42	18	70	1	S113011
50	51	22	90	1	B103008	50	51	18	85	1	S113002
100	64	18	105	1	B103009	100	64	18	105	1	S113006
150	74	28	118	1		150	74	28	118	1	S113010
200	79	34	135	1	B103011	200	79	34	131	1	S113012
250	85	34	145	1	B103012	250	85	34	140	1	S113014
300	87	34	160	1	B103013	300	87	34	156	1	S113016
500	105	34	180	1	B103016	500	105	34	180	1	S113020
1 000	131	42	220	1	B103020	1 000	131	36	220	1	S113026
2 000	166	50	280	1	B103022	2 000	166	45	280	1	S113030
3 000	187	50	310	1	B103023	3 000	187	46	310	1	S113032
5 000	220	50	365	1	B103025	5 000	220	46	365	1	S113036