

FAD Na₂

flavin-adenin-dinukleotid disodná sůl, riboflavin-5'-adenosin difosfát disodná sůl


CAS: 84366-81-4

Vzorec: C₂₇H₃₁O₁₅N₉Na₂P₂

Technické parametry

Molekulová hmotnost:

865,55



Čistota	Balení	Katalogové číslo
min 95 %, pro biochemii	10 mg	R.5581.2
min 95 %, pro biochemii	25 mg	R.5581.3
min 95 %, pro biochemii	100 mg	R.5581.1
min 95 %, pro biochemii	500 mg	R.5581.4
min 95 %, pro biochemii	1 g	R.5581.5
min 95 %, CELLPURE®	10 mg	R.6833.1
min 95 %, CELLPURE®	50 mg	R.6833.2
min 95 %, CELLPURE®	100 mg	R.6833.3