

Kádinka nízká s výlevkou

Kádinka z borosilikátového skla BORO 3.3 nebo skla PYREX®



Kádinka nízká s výlevkou PYREX



Kádinka nízká s výlevkou BORO 3.3

BORO 3.3					PYREX®				
Objem (ml)	Průměr (mm)	Výška (mm)	Balení (ks)	Katalogové číslo	Objem (ml)	Průměr (mm)	Výška (mm)	Balení (ks)	Katalogové číslo
5	22	30	1	T100703	5	22	30	1	S100003
10	26	35	1	T100704	10	26	35	1	S100005
25	34	50	1	T100706	25	34	50	1	S100013
50	42	60	1	T100708	50	42	60	1	S100002
100	50	70	1	T100709	100	50	70	1	S100004
150	60	80	1	T100710	150	60	80	1	S100006
250	70	95	1	B100712	250	70	95	1	S100010
400	80	110	1	B100715	400	80	110	1	S100014
600	90	125	1	B100717	600	90	125	1	S100018
800	100	135	1	B100719	800	100	135	1	S100023
1 000	105	145	1	B100720	1 000	105	145	1	S100024
2 000	130	185	1	B100722	2 000	130	185	1	S100028
3 000	150	210	1	B100723	3 000	150	210	1	S100030
5 000	170	270	1	B100225	5 000	170	270	1	S100034
10 000	217	350	1		10 000	217	350	1	S100038