

Nuč Filtrační

Filtrační nuč z borosilikátového skla

- Autoklárovatelná

Technické parametry

Materiál:

borosilikátové sklo



Nuč filtrační

Objem (ml)	Porozita	Průměr nálevky (mm)	Průměr stonku (mm)	Balení (ks)	Katalogové číslo
50	1	35	10	1	R.Y772.1
50	2	35	10	1	R.Y773.1
50	3	35	10	1	R.Y774.1
50	4	35	10	1	R.Y775.1
50	5	35	10	1	R.Y871.1
75	1	45	10	1	R.Y777.1
75	2	45	10	1	R.Y778.1
75	3	45	10	1	R.Y779.1
75	4	45	10	1	R.Y780.1
75	5	45	10	1	R.Y872.1
125	1	60	10	1	R.Y782.1
125	2	60	10	1	R.Y783.1
125	3	60	10	1	R.Y784.1
125	4	60	10	1	R.Y785.1
125	5	60	10	1	R.Y873.1
500	1	95	22	1	R.Y787.1
500	2	95	22	1	R.Y788.1
500	3	95	22	1	R.Y789.1
500	4	95	22	1	R.Y790.1
500	5	95	22	1	R.Y874.1

1000	1	120	22	1	R.Y792.1
1000	2	120	22	1	R.Y793.1
1000	3	120	22	1	R.Y794.1
1000	4	120	22	1	R.Y795.1
1000	5	120	22	1	R.Y875.1
4000	1	175	30	1	R.Y876.1
4000	2	175	30	1	R.Y877.1
4000	3	175	30	1	R.Y878.1
4000	4	175	30	1	R.Y879.1