

Nuč s rovnou fritou | SIMAX®

Filtrační nuč s fritou různé pórovitosti

Technické parametry

Porozita	Velikost pórů (μm)	Značení dle ISO 4793
S1	100 - 160	P 160
S2	40 - 100	P 100
S3	16 - 40	P 40
S4	10 - 16	P 16



Objem (ml)	Porozita	Průměr frity (mm)	Výška nálevky (mm)	Balení (ks)	Katalogové číslo
35	S1	30	143	1	S334103
35	S2	30	143	1	S334203
35	S3	30	143	1	S334303
35	S4	30	143	1	S334403
70	S1	40	165	1	S334104
70	S2	40	165	1	S334204
70	S3	40	165	1	S334304
70	S4	40	165	1	S334404
150	S1	60	191	1	S334106
150	S2	60	191	1	S334206
150	S3	60	191	1	S334306
150	S4	60	191	1	S334406
500	S1	90	246	1	S334109
500	S2	90	246	1	S334209
500	S3	90	246	1	S334309
500	S4	90	246	1	S334409
1000	S1	120	306	1	S334112
1000	S2	120	306	1	S334212

1000	S3	120	306	1	S334312
1000	S4	120	306	1	S334412
2500	S1	150	426	1	S334115
2500	S2	150	426	1	S334215
2500	S3	150	426	1	S334315
2500	S4	150	426	1	S334415