

Podnos - miska ABS | KARTELL

Tác pro laboratorní použití

Technické parametry

Materiál:

bílý ABS



Podnos - miska ABS

Celkové rozměry (mm)	Rozměry dna (mm)	Hloubka (mm)	Balení (ks)	Katalogové číslo
201 x 151 x 21	170 x 110	10	1	K005701
303 x 151 x 21	270 x 120	10	1	K005702
353 x 252 x 21	320 x 220	10	1	K005703
254 x 254 x 21	220 x 220	10	1	K005704
408 x 300 x 21	380 x 270	10	1	K005713
201 x 151 x 41	150 x 100	30	1	K005705
252 x 252 x 42	210 x 210	30	1	K005708
303 x 151 x 42	260 x 110	30	1	K005706
353 x 253 x 41	310 x 210	40	1	K005707
409 x 300 x 42	360 x 250	40	1	K005714
202 x 151 x 81	140 x 90	70	1	K005709
252 x 252 x 82	200 x 200	70	1	K005712
303 x 151 x 81	240 x 90	70	1	K005710
352 x 252 x 81	290 x 190	70	1	K005711
408 x 299 x 81	340 x 240	70	1	K005715